EN Do not reuse!
ES No reutilizar!
PT Não reutilizar!
IT Non riutilizzare!
FR Ne pas réutiliser!
DE Nicht reesterilisieren!
NL Niet opnieuw steriliseren!
TR Yeniden sterilize etmeyin!
EL Δεν επαναποστειρώνεται!
PL Nie nadaje się do powtórzonego stosowania!
DA Må ikke genbruges!
SR Samo za jednokratno korišćenje!
CZ Nepoužívat opakovaně!
SK Nepoužívať opakovaně!
HR Samo za jednokratnu uporabu!

EN Store between 5 °C and 25 °C
ES Conservar a una temperatura de 5 °C a 25 °C
PT Armazenar entre 5 °C e 25 °C
IT Conservare ad una temperatura fra 5 °C e 25 °C
FR Conserver a une température comprise entre 5 °C et 25 °C
DE Bei 5 °C - 25 °C lagern
NL Tussen 5 °C en 25 °C bewaren
TR 5 °C ve 25 °C arasında saklayın
EL Θυλάξτε σε θερμοκρασία μεταξύ 5°C και 25°C
PL Przechowywać w temperaturze 5 °C – 25 °C
DA Opbevares ved temperaturer mellem 5 °C og 25 °C
SR Skladati se između 5 °C i 25 °C
CZ Uchovávejte při teplotě 5 °C až 25 °C
SK Skladujte medzi 5 °C a 25 °C
HR Skladati između 5 °C i 25 °C

EN Expiry Date, Year/Month
ES Fecha de caducidad, año/mes
PT Data de expiração, ano/mes
IT Data di Scadenza, Anno/Mese
FR Date d'expiration, annee/mois
DE Haltbarkeitsdatum (Jahr/Monat)
NL Houdbaarheidsdatum (Jaar/Maand)
TR Son Kullanma Tarihi, Yıl/Ay
EL Ημερομηνία λήξης - Έτος/Μηνός
PL Data Ważności, Rok/Miesiąc
DA Udløbsdato, år/måned
SR Datum isteka roka trajanja, godina/mjesec
CZ Datum ukončení doby použitelnosti, rok/měsíc
SK Datum expirácie, rok/mesiac
HR Datum isteka roka trajanja, godina/mjesec

EN Do not use if packaging is damaged
ES No utilizar si el envase está dañado
PT Não usar de a embalagem estiver danificada
IT Non utilizzare se la confezione è danneggiata
FR Ne pas utiliser si l'emballage est endommagé
DE Bei beschädigter Verpackung nicht verwenden
NL Niet gebruiken als de verpakking beschadigd is
TR Ambalaj hasarlıysa kullanmayın
EL Μην χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά
PL Nie używać, gdy opakowanie jest uszkodzone
DA Må ikke anvendes, hvis emballagen er beskadiget
SR Ne koristiti ako je pakovanje oštećeno
CZ Nepoužívejte, pokud je obal poškozen
SK Nepoužívať, ak je obal poškodený
HR Ne upotrebljavati ako je ambalaza oštećena

EN Do not re-sterilize
ES No reesterilizar
PT Não re-esterilizar
IT Non riesterilizzare
FR Ne pas reesteriliser
DE Nicht reesterilisieren
NL Niet opnieuw steriliseren
TR Yeniden sterilize etmeyin
EL Δεν επαναποστειρώνεται
PL Nie steryliżować ponownie
DA Må ikke reesteriliseres
SR Ne sterilizovati više puta
CZ Opakované nesterilizuje
SK Nepopakovať sterilizáciu
HR Ne sterilizirati više puta

EN Date of Manufacture, Year/Month
ES Fecha de fabricación, año/mes
PT Data de fabricação, ano/mes
IT Data di fabbricazione, Anno/Mese
FR Date de fabrication, annee/mois
DE Herstellungsdatum (Jahr/Monat)
NL Productiedatum (Jaar/Maand)
TR Üretim Tarihi, Yıl/Ay
EL Ημερομηνία κατασκευής - Έτος/Μηνός
PL Data Produkcji, Rok/Miesiąc
DA Produktionsdato, år/måned
SR Datum proizvodnje, godina/mjesec
CZ Datum výroby, rok/měsíc
SK Datum výroby, rok/mesiac
HR Datum proizvodnje, godina/mjesec

EN Protect from humidity
ES Proteger de la humedad
PT Proteger da umidade
IT Proteggere dall'umidità
FR Conserver a l'abri de l'humidité
DE Vor Feuchtigkeit schützen
NL Droog bewaren
TR Nemden koruyun
EL Να αποφεύγετε την υγρασία
PL Chronić przed wilgocią
DA Beskyttes mod fugt
SR Zaštiti od vlage
CZ Chráněte před vlhkem
SK Chránite pred vlhkom
HR Zaštiti od vlage

EN Batch Number
ES Código de lote
PT Lote nº
IT Numero di Lotti
FR Numéro de lot
DE Chargennummer
NL Lotnummer
TR Parti Numarası
EL Αριθμός Παρτίδας
PL Numer Partii
DA Batchnummer
SR Broj šarže
CZ Číslo šarže
SK Kód šaržky
HR Broj šarže

EN Read IFU before use
ES Leer las instrucciones antes de usar
PT Ler as instruções antes do uso
IT Leggere le istruzioni prima dell'uso
FR Lire les instructions d'utilisation avant toute utilisation
DE Vor Verwendung Gebrauchsanleitung lesen
NL Lees de gebruiksaanwijzing voor gebruik
TR Kullanmadan önce KTY'i okuyun
EL Ανοίξτε τις οδηγίες χρήσης πριν το χρησιμοποιήσετε
PL Przed użyciem zapoznać się ze Wskazówkami dotyczącymi stosowania produktu
DA Må ikke reesteriliseres
SR Ne sterilizovati više puta
CZ Opakované nesterilizuje
SK Nepopakovať sterilizáciu
HR Ne sterilizirati više puta

EN Keep away from sunlight
ES Mantener alejado de la luz del sol
PT Manter longe da luz do sol
IT Tenere lontano dalla luce diretta
FR Conserver a l'abri de la lumière du soleil
DE Vor Sonneneinstrahlung schützen
NL Uit de zon bewaren
TR Güneş ışığından uzak tutun
EL Να μην εκτίθεται στον ήλιο
PL Chronić przed działaniem promieni słonecznych
DA Må ikke udsættes for direkte sollys
SR Čuvati zaštićeno od sunčeve svetlosti
CZ Chráněte před sluncem
SK Chránite pred slnkom
HR Čuvati zaštićeno od sunčeve svetlosti

EN Reference, Product code
ES Referencia, código del producto
PT Referência, Código do Produto
IT Riferimento, Codice del Prodotto
FR Référence, code du produit
DE Referenz, Produktcode
NL Referentie, productcode
TR Referans, Ürün kodu
EL Αναφορά, κωδικός προϊόντος
PL Numer referencyjny, Kod produktu
DA Referencia, produktkode
SR Referenca, kod proizvoda
CZ Reference, kód výrobku
SK Referencia, katalógové číslo
HR Referencija, šifra proizvoda

EN Sterile R: Gamma Irradiation
ES Esterilizado mediante irradiación gamma
PT Estéril R: irradiação gama
IT Sterile R: irradiazione Gamma
FR Sterile R: irradiation aux rayons gamma
DE Steril R: Gammastrahlen
NL Steriel R: Gammastralen
TR Steril R: Gama Iyrtılması
EL Αποστειρωμένο
PL Produkt jałowy: sterylizowany promieniowaniem gamma
DA Steril R: Gammabestråling
SR Sterilno R: gama zračenje
CZ Sterilizováno gama zářením
SK Sterilizované žiarením
HR Sterilno R: gama zračenje

INSTRUCTIONS FOR USE CuraMedical B.V. - We help to cure people

CuraCel® Sterile absorbable haemostat, type Oxidised Regenerated Cellulose (ORC)

CuraCel® Standard

CuraCel® High Density

CuraCel® Fibrillar

- EN** Instructions for Use
- ES** Instrucciones de Uso
- PT** Instruções de utilização
- IT** Istruzioni per l'uso
- FR** Instructions d'utilisation
- DE** Gebrauchsanleitung
- NL** Gebruiksaanwijzing
- TR** Kullanma Talimatı
- EL** ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- PL** Wskazówki dotyczące stosowania opatrunku
- DA** Brugsvejledning
- SR** Uputstvo za korišćenje
- CZ** Navod k použití
- SK** Navod na použitie
- HR** Uputa za uporabu

CuraMedical B.V.
 Industrieweg 6 B, 1566 JP Assendelft
 Industrieweg 13-1, 1566 JN Assendelft
 The Netherlands
 Tel: +31 75 657 1885 www.curamedical.com
 Fax: +31 75 657 1184 sales@curamedical.com
 IFU-CuraCel-REV06a-2017-12

CuraCel® Standard ABSORBABLE HAEMOSTAT (oxidised regenerated cellulose)

Labelling:
Medical Cover

Labelling :
Individual (Sachet)
Pack / Sterile Pack

Labelling :
Box (Protective Cover) - Front

CuraCel® High Density ABSORBABLE HAEMOSTAT (oxidised regenerated cellulose)

Labelling:
Tyvek Cover

Labelling :
Individual (Sachet)
Pack / Sterile Pack

Labelling :
Box (Protective Cover) - Front

CuraCel® Fibrillar ABSORBABLE HAEMOSTAT (oxidised regenerated cellulose)

Labelling:
Tyvek Cover

Labelling :
Individual (Sachet)
Pack / Sterile Pack

Labelling :
Box (Protective Cover) - Front

Index

EN Instructions for Use	1
ES Instrucciones de Uso	2
PT Instrucoes de utilizacao	3
IT Istruzioni per l'uso	4
FR Instructions d'utilisation	5
DE Gebrauchsanleitung	6
NL Gebruiksaanwijzing	7
TR Kullanma Talimati	8
EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	9
PL Wskazowki dotyczące stosowania opatrunku	10
DA Brugsvejledning	11
SR Uputstvo za korišćenje	12
CZ Navod k použití	13
SK Navod na použitie	14
HR Uputa za uporabu	15

CuraCel® Standard

CuraCel® High Density

CuraCel® Fibrillar

EN INSTRUCTIONS FOR USE

ABSORBABLE HAEMOSTAT (oxidised regenerated cellulose)

DESCRIPTION

Curacel® Absorbable Haemostat is a sterile product made of oxidised regenerated cellulose (polyoxyanhydro glucuronic acid) available in the flat form with textile structure or in fibrillose form. Curacel® Haemostat can be cut into the desired dimension. Its structure is stable and should be stored at controlled room temperature. The fabric is pale yellow and has a faint, caramel-like aroma. A slight discoloration may occur with time, but this does not affect its performance.

ACTIONS

Curacel® Haemostat is applied dry. It can be laid on or held against a bleeding surface. In vascular surgery, Curacel® Haemostat can be very loosely wrapped around a bleeding surface (see precautions and adverse reactions). As special features, Curacel High Density can be sutured, and Curacel Fibrillar can be peeled off in layers for optimal dosage. Action mechanism of Curacel® is independent from blood coagulation mechanism of the body. When it contacts with blood, Curacel® Haemostat forms a brownish or black gelatinous mass which aids clot formation in 3-4 minutes. This gelatinous mass acts as a physical matrix to which platelets can adhere. With platelet aggregation and formation of platelet-fibrin plug, haemostasis occurs. When Curacel® Haemostat is used properly in minimal amounts, it is absorbed from the sites of implantation without tissue reaction. Absorption depends on; Haemostat amount used, tissue type and blood saturation degree. Curacel® should be saturated completely with blood to obtain the absorption effect. It is completely absorbed within 7-14 days.

INDICATIONS

Curacel® Haemostat is designed to arrest capillary bleeding and bleeding from parenchymatous organs and resection areas at surgical interventions. It shall be used to stop bleeding when other haemostatic procedures are ineffective or impractical. It is suitable for use in general surgery and digestive surgery, neurosurgery (especially cerebral operations), plastic surgery, orthopaedic, gynaecology, urology, stomatology, traumatology, and many other branches of surgery. Curacel® Haemostat can be applied into cavities (after extirpation of tumours) as well as endoscopic interventions or dental praxis.

CONTRAINDICATIONS

Curacel® Haemostat should not be used for packing. If you have to use in this manner, it must be removed after haemostasis is achieved. Curacel® Haemostat should not be used for implantation in bone defects, such as fractures, since it may interfere with callus formation and may cause cyst formation. Curacel® Haemostat should not be used to control haemorrhage from large arteries. Since body fluids, except whole blood, such as serum, don't react with Curacel® Haemostat, it should not be used to produce haemostatic effect on non-haemorrhagic serous oozing surface.

WARNINGS

Curacel® Haemostat should not be re-sterilized, not even by gamma-irradiation, because product conformity cannot be guaranteed. Closing Curacel® Haemostat in a contaminated wound without drainage may lead to complications and should be avoided. The haemostatic effects of Curacel® Haemostat are greater when it is applied dry; it should not be moistened with water or saline prior to use. Impregnation of Curacel® Haemostat with buffering materials or Haemostatic substances should not be done. Curacel® Haemostat should not be used together with anti-infective agents including antibiotics. Since thrombin activity is destroyed due to the low pH of the product, haemostatic effect of Curacel® Haemostat is not enhanced by addition of thrombin. Curacel® Haemostat may be left in situ when necessary, but it is advisable to remove it once haemostasis is achieved. It must always be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm regardless of the type of surgical procedure, because Curacel® Haemostat may exert pressure resulting in paralysis and/or nerve damage by swelling. Re-motion of Curacel® Haemostat could possibly occur by means such as repacking, further intraoperative manipulation, lavage, etc. There have been reports that Oxidised (Regenerated) Cellulose Haemostat, when left in the patient after closure during procedures such as lobectomy, laminectomy and repair of a frontal skull fracture and lacerated lobe, migrated from the site of application into foramina in bone around the spinal cord resulting in paralysis and, in another case, migrated to the left orbit of the eye causing blindness. While these reports cannot be confirmed, physicians must act cautiously, regardless of the type of surgical procedure, and consider the advisability of removing Curacel® Haemostat once haemostasis is achieved. Curacel® Haemostat should not be used as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections. Curacel® Haemostat should not be used as an adhesion prevention product.

PRECAUTIONS

Use only minimal amounts of Curacel® Haemostat in appropriate size needed to obtain haemostasis, holding it firmly in place until bleeding stops. Remove any excess before closure in order to facilitate absorption and minimize the possibility of foreign body reaction. In urological procedures, minimal amounts of Curacel® Haemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product. Since absorption of Curacel® Haemostat could be prevented in chemically cauterized areas, its use should not be preceded by application of silver nitrate or any other escharotic chemicals. If Curacel® Haemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges. It should also be removed from

open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped. Precautions should be taken in otorhinolaryngology surgery to assure that none of the material is aspirated by the patient. (Examples: controlling haemorrhage after tonsillectomy and controlling epistaxis). Care should be taken not to apply Curacel® Haemostat tightly when it is used as a wrap during vascular surgery.

ADVERSE REACTIONS

"Encapsulation" of fluid and foreign body reactions has been reported for Oxidised (Regenerated) Cellulose Haemostat. Infection post neck surgery has been reported due to misuse of Oxidised (Regenerated) Cellulose (too much material used, no complete blood saturation and the product was layered). It is important to use the product in conformance with the IFU. There have been reports of stenotic effect when Oxidised (Regenerated) Cellulose Haemostat has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of Oxidised (Regenerated) Cellulose Haemostat, it is important to be cautious and avoid applying the material tightly as a wrapping. Paralysis and nerve damage have been reported when Oxidised (Regenerated) Cellulose Haemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm. While most of these reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures. Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when Oxidised (Regenerated) Cellulose Haemostat was placed in the anterior cranial fossa. Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy has been reported. There has been one report of a blocked ureter after kidney resection, in which post-operative catheterization was required. Occasional reports of "burning" and "stinging" sensation and sneezing when Oxidised (Regenerated) Cellulose Haemostat has been used as packing in epistaxis, are believed to be due to the low pH of the product. Burning has been reported when Oxidised (Regenerated) Cellulose Haemostat was applied after nasal polyp removal and after haemorrhoidectomy. Headache, burning, stinging, and sneezing in epistaxis control and other rhinological procedures, and stinging when Oxidised (Regenerated) Cellulose Haemostat was applied on surface wounds (varicose ulcerations, dermabrasions, and donor sites) have been reported as well.

DOSAGE AND ADMINISTRATION

Sterile technique should be observed in removing Curacel® Haemostat from its sterile container. Minimal amount of Curacel® Haemostat in appropriate size is laid on the bleeding site or held firmly against the tissues until haemostasis is achieved. Required amount of Curacel® Haemostat depends on the nature and intensity of the haemorrhage to be stopped. The haemostatic effect of Curacel® Haemostat is optimal when used dry. Moistening the material with water or physiological saline solution prior to use is not recommended.

CAUTION

Since Curacel® Haemostat cannot be re-sterilized, opened and unused Curacel® Haemostat should be discarded. Contaminated and/or excess material must be disposed in accordance with the applicable hospital guidelines for the disposal of clinical and related waste.

NOTES AND INFORMATION ON SHELF-LIFE

Curacel® Haemostat should be stored dry under controlled conditions (5°C – 25°C) in the original packaging, and protected from direct sunlight. Shelf life of the product is 3 years. The expiry date of Curacel® Haemostat is printed on the pack. Do not use Curacel® Haemostat after this date.

PRESENTATION AND CONTENTS

Curacel® Haemostat is available sterile, as an original pack in different sizes. Curacel® Haemostat strips are sterile-packed individually.

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