# CuraSpon<sup>®</sup>

LCS003.07

<b>EN</b>	Instructions for use	4
<b>NL</b>	Gebruiksaanwijzing	6
<b>DA</b>	Brugsvejledning	8
<b>DE</b>	Gebrauchsinformation	10
<b>FR</b>	Notice d'utilisation	12
<b>ES</b>	Instrucciones de uso	14
<b>PT</b>	Instrucoes para utilizacco	16
<b>TR</b>	Kullunim talimati	18
<b>PL</b>	Instrukcje użytowania	20
<b>IT</b>	Spugna di gelatina riassorbibile	22
<b>EL</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	24
<b>CZ</b>	Návod k použití houby	26
<b>SK</b>	Návod na použitie špongie	28
<b>RO</b>	Instrucțiuni pentru utilizare	30
<b>SR</b>	Uputstvo za korišćenje	32
<b>HR</b>	Uputa za uporabu	34

## 1. Product Description

The CuraSpon® gelatine sponge is made of 100% purified porcine gelatine. It is intended for the use as an absorbable haemostat in various fields of surgery. It is insoluble in water and fully absorbable. Its porous structure and interstices enables the product to absorb more than 35 times its own weight of blood and fluids. The absorptive capacity of CuraSpon is a function of its physical size, increasing as the amount of gelatine increases. In addition, because of its neutral character and nature, it may be used as a carrier for any drug (the product literature for the applicable drug should be consulted for complete prescribing information). It is supplied in several sizes. It is packaged individually in (double) blister, pouch or jar/container packaging, sterilized by gamma rays, and is for single use only.

## 2. Indications

- 2.1 General indications: the CuraSpon® gelatine sponge is used during and after surgical procedures to achieve local haemostasis by tamponade effect.
  - 2.2 In Dental surgery, the CuraSpon® cube is an aid in providing haemostasis and obliterating dead space created by extraction of the teeth, root amputations and removal of cysts, tumours and impacted teeth.
  - 2.3 For ENT surgery the use has been proven to be effective in supporting and keeping the fascia or perichondrium in the middle ear (dry application) in place or in the outer ear canal after tympanoplasty.
  - 2.4 For rectal surgery, haemorrhoid operations and gynaecology CuraSpon® gelatine sponge is used to control bleeding.
- ## 3. Contraindications
- Allergies to porcine products.
- The gelatine sponge should not be used in closures of skin incisions, as it may interfere with wound healing.

- The product should not be used without antibiotics in infected wounds.

## 4. Side effects

There have been no reported adverse effects for product used correctly according to these Instructions For Use and when over-packing is not applied. Formation of tissue granulation during ENT/middle ear procedures has been reported in at least one animal study.

## 5. Directions for use:

The gelatine sponge can be used dry or saturated with a sterile physiological saline solution.

If used dry, the sponge is cut into the desired size and is slightly compressed. The sponges must be applied to a bleeding area under light pressure for one or two minutes until the bleeding stops.

When used with saline, CuraSpon should be soaked in the solution, then withdrawn, squeezed thoroughly between gloved fingers to expel air bubbles present in the interstices, replaced in saline, and kept there until needed.

If bleeding is controlled after (possible multiple) applications of the gelatine sponge, the material can be left in situ. The excess should be removed.

Depending on the method of use and the operated site, the total absorption time of the sponge varies from four to six weeks. When applied to bleeding nasal, rectal or vaginal mucosa, it liquefies within 2 to 5 days. It should be noted that the absorption is a process where the volume as well as the absorption characteristics continuously decrease over time.

- 5.1 Gelatine Powder (powder paste); the powder is prepared by milling absorbable gelatine sponges. One (1) gram of the powder is mixed with 3 to 4 ml of a sterile liquid solution, then kneaded with the gloved fingers until the desired level of viscosity and consistency is reached. It is indicated in procedures involving cancellous bone bleeding. The paste may be smeared or pressed against the bleeding surface to control bleeding. Accurate application of the paste to the bleeding site can be achieved by using a syringe. When bleeding stops the excess should be removed.

- 5.2 ENT applications: The gelatine sponge is cut as needed to fit the otic, nasal or oral cavity or defect and inserted to support and separate tissues to prevent adhesions and/or to control bleeding by tamponade effect. The excess should be removed.

## Important Notice!

In cases where the gelatine sponge has been compressed (either by the manufacturer or end-user), the gelatine sponge may expand with the absorption of fluids (the sponge will obtain its original shape). Where possible CuraSpon® should be removed after usage and bleeding has stopped in radical cavities, laminectomy procedures, around or in proximity to foramina in bone, areas of bony confine, the spinal cord and/or the optic nerve and chiasm or closed tissue spaces with presence of bone. This might lead to unintended pressure on neighboring structures which may result in pain for the patient or might create the potential for nerve damage.

The over-packing of CuraSpon®, particularly within bony cavities, should be avoided, since recovering to its initial volume may interfere with normal function and/or could cause possible or eventual compression necrosis of surrounding tissue. Formation of tissue granulation has been observed during ENT / middle ear procedures.

In cases of postoperative infection, re-operation may be necessary to remove infected material and allow drainage.

## Warnings

Do not use when the package has been opened or has been damaged.

The absorbable gelatine CuraSpon® sponge is a single-use product which is not suitable for re-sterilization.

## 6. General storage and handling

The product should be stored in its sealed protective outer packaging. Prior to opening, the packaging ensuring the integrity of the sterile barrier has to be inspected. Open the package just prior to use in a patient.

- 6.1 Packaging and sterility  
The products are packaged sterile.  
The packaging consists of:  
- Storage and transport packaging (Secondary packaging)  
- Sterile packaging (primary packaging)  
The packaging is subject to and conforms to European & International legislation and other applicable standards. The packaging protects the product from external influences and guarantees its sterility during storage.
- 6.2 Handling of the sterile packaging  
Remove the product from the sterile package using acceptable aseptic technique.
- 6.3 Sterilization  
The products are sterilized by Gamma radiation.
- 6.4 Storage  
The products should be stored in a dry place in its protective packaging at room temperatures. Avoid contact with direct sunlight.

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